

## Tools & Supplies

### Tools


- 2 x #5 Forceps Fine-tip metal forceps
- 2 x #2 Forceps Medium-tip metal forceps
- 1 x #7 Forceps Curved medium-tip metal forceps
- 1 x Non-Magnetic Epoxy-Coated AA Forceps Medium forceps, #2AP
- 1 x Surgical Scissors 10.5 cm, straight blade
- 1 x Vannas scissors, 3" Small dissecting spring scissors, straight
- 1 x Scalpel Handle
- 1 x Retraction Kit for small animals from FST
- 1 x Hemostats 6" Straight or curved For inserting / removing Omnetics connectors

### Supplies

- 6-0 Sutures, Nylon monofilament; Used to close incision site
- 8-0 Sutures, Nylon monofilament; Used to secure array to muscle
- Dental Cement; Used to adhere Omnetics connector to skull
- Sterile Saline; Used to keep muscle pliable
- Scalpel Blades, #10 or #11; Used to score skull for improved adhesion

## Hardware

Myomatrix arrays are equipped with a high-density connector:



The arrays are assembled with a 36-pin Omnetics connector with 4 guide posts — Omnetics part number NPD-36-VV-GS (A79026) which is compatible with a variety of electrophysiology recording hardware. Note that CAMBER recommends for users new to EMG to record using a **bipolar headstage**. Bipolar mode for recording can give better stability for the recordings.

Below is a list of headstages that are compatible with the Omnetics NPD-36-VV-GS connector based on the descriptions available online, though not an exclusive list. Compatibility with the 36-pin Omnetics connector is the key consideration in choosing an amplifying system.

Last edited  
04/29/2024

| <b>Company</b>           | <b>Option 1</b>          | <b>Option 2</b>               | <b>Option 3</b> | <b>Option 4</b>      |
|--------------------------|--------------------------|-------------------------------|-----------------|----------------------|
| <b>Blackrock</b>         | CerePlex mu (32 channel) |                               |                 |                      |
| <b>Intan, Open Ephys</b> | RHD (32 or 64 channels)  | RHS 32                        |                 |                      |
| <b>NeuraLynx</b>         | HS-64-mux-PTB            | QuickClip 72 to Dual Omnetics | HS-36-PTB       | ADPT-HS36-N2<br>T-32 |
| <b>Ripple</b>            | Nano2                    | Pico (32 channel)             |                 |                      |
| <b>Plexon</b>            | F3                       | R                             |                 |                      |